

## Uvex Superfit Clean Room Safety Spectacle - Pack of 5



### Why Choose Uvex Superfit Clean Room Safety Spectacle - Pack of 5

Lightweight, sterile safety eyewear weighing 21g  
Suitable for cleanroom environments  
Extremely good coverage of brow and cheek area  
Comfortable, non-slip fit  
First autoclavable protective spec with anti-fog lens coating  
Autoclavable at least 10 times  
Pack size: 5 (=1 x unit)  
Uvex code: 9178500

### More Info

Lightweight safety spec weighing 21g  
Suitable for cleanroom environments  
Extremely good coverage of brow and cheek area  
Comfortable, non-slip fit  
First autoclavable protective spec with anti-fog coating  
XST side arms for increased comfort and pressure free fit  
Wide panoramic field of vision  
Professional protection of products and humans for medical and pharmaceutical industries  
Sterile eye protection in every situation  
Autoclaveable at least 10 times at 121°C for 20 minutes  
Innovative design for maximum functionality and wearer comfort  
Non-allergenic material which is free of harmful substances  
Ultra-thin wraparound lenses  
Applications: Medical and pharmaceutical industries, laboratories, clinics, cleanroom environments, food industries

Disinfection: Commercially available disinfectants can be used, either by spraying or wiping. Goggles should not be placed in disinfectant for an extended period of time.

## Brand

Uvex

## Type

Goggles/Specs

## Lead time (if out of stock)

10-15 days

## Pack Size

5

## Size

One size

## Colour

Clear

## Applications

Medical and pharmaceutical industries, laboratories, clinics, cleanroom environments, food industries

## Technical data

Frame front piece/temple: PC blend, white  
Frame temple end: PU, white/light blue  
Frame EN Norm: DIN166  
Lens material: PC  
Lens thickness: 1.5mm  
Lens coating: supravision CR (anti-fog on the inside, scratch resistant on the outside)  
Lens colour: clear  
Lens filter: 2C-1,2  
Lens optical class: 1  
Lens impact resistance: F (45 m/s)  
Lens surface resistant to fine particles: K  
Lens antifog: N

## Packaging

5 goggles per box

## Sterility

Sterile

## CE Certified

Yes

## Disposal

After 10 autoclave cycles, the spectacle should be replaced

## Washing

Commercially available disinfectants can be used, either by spraying or wiping. Goggles should not be placed in

disinfectant for an extended period of time.



[www.cleanroomshop.com](http://www.cleanroomshop.com)  
Tel: + 44 (0) 1524 813022

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)  
Fax: + 44 (0) 1524 811589