

Why Choose Disposable Latex Gloves 16" Sterile - Extra

BioClean Gloves - Extra

Latex & textured grip - extremely tactile

Individually double bagged & sterile

ISO Class 4+ compatible

Length: 400mm (16")

Pack size: 10 pairs

Colour: natural

Full validation pack available

Nitritex code: BLAS

Available Options

More Info

At 400mm (16") BioClean Extra™ Sterile Latex Gloves provide exceptional comfort and elbow length protection. With a textured surface for enhanced grip, this natural-coloured glove is ambidextrous and features a beaded cuff for stability on the arm.

Disposable Latex Gloves 16" Sterile - BioClean - Extra

Length: 400mm (16") length

Offer elbow length protection and are powder free

Each pair is double bagged in 'Easy Tear' plastic packaging that is silicone, amine and amide-free

Gamma sterilised and suitable for use in sterile ISO class 4 cleanroom environments

The beaded cuff ensures a secure fit at the wrist

Textured surface provides a good grip even in damp or wet operating environments

Antistatic: Not static dissipative but Latex is an excellent insulator



www.cleanroomshop.com

info@cleanroomshop.com

Tel: + 44 (0) 1524 813022

Fax: + 44 (0) 1524 811589

Large

XL

Product ordering

Size: Small (GLTS2S), Medium (GLTS2M), Large (GLTS2L), XL (GLTS2XL)

Product Icon Attributes

Suitable for ISO 4-8, Individually wrapped , 400mm length, 10 pairs

Brand

BioClean

Lead time (if out of stock)

2-3 weeks

Pack Size

10 hand-specific pairs

Construction/Material

Latex

Length

400mm (16")

Palm Size

Palm widths (±5mm)

Small: 85mm

Medium: 95mm

Large: 106mm

X Large: 115mm

Test method: EN455-2:2009

Colour

Natural

Applications

ISO Class 4+ compatible. The gloves have been developed for users who require increased tactility and control against product contamination.

Technical data


Complies with the requirements of the Glove Standard EN420:2003, the European Medical Gloves Standard EN455 (Parts 1, 2 & 3) & EN374-1, 2 & 3:2003. Manufactured in a facility holding the ISO 9001:2008 and 89/686/EEC Article 11b certifications. Processed in an NEBB certified ISO Class 5 cleanroom environment.

Packaging

Double bagged in 'Easy Tear' plastic packaging that is silicone, amine and amide-free

Sterility

Sterile



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589

Handling & storage

Store in a dry cool place (<40°C) away from direct sunlight and fluorescent light

Disposal

Recommended for single use



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589