



Why Choose Disposable Latex Gloves 16" Non Sterile - Legion

BioClean Gloves - Legion

Latex & textured grip - extremely tactile

Double bagged & cleanroom packed

ISO Class 4+ compatible

Length: 400mm (16")

Pack size: 100 ambidextrous gloves

Colour: natural

Full validation pack available

Nitritex code: BLA-3

Available Options



More Info

For supreme comfort, exceptional flexibility and elbow length protection, the 400mm (16") BioClean Legion™ Latex Glove is the natural choice. With a textured surface for enhanced grip, the ambidextrous natural-colour glove is powder-free.

Disposable Latex Gloves 16" Non Sterile - BioClean - Legion BLA3

Length: 400mm (16") elbow length protection

Powder free, ISO Class 4 compatible and double bagged

The textured surface provides a good grip even in damp or wet operating environments

Latex provides exceptional flexibility and comfort

Contains a latex protein level of less than 50µg/gram

Textured for improved grip

Processed to ensure ISO Class 4 compatibility

Beaded cuff provides easier donning and stability on the arm

Antistatic: Not static dissipative but Latex is an excellent

www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

sales@cleanroomshop.com
Fax: + 44 (0) 1524 813022

Medium

Large

XL

Product ordering

Sizes: Small (GLTNS2S), Medium (GLTNS2M), Large (GLTNS2L), XL (GLTNS2XL)

Product Icon Attributes

Suitable for ISO 4-8, Double bagged pack, 100 gloves, 400mm length

Brand

BioClean

Lead time (if out of stock)

2-3 weeks

Pack Size

100 ambidextrous gloves

Construction/Material

Latex

Length

400mm (16")

Colour

Natural

Applications

ISO Class 4+ compatible. The gloves have been developed for users who require increased protection and control against product contamination.

Technical data

Complies with the requirements of the Glove Standard EN420:2003, the European Medical Gloves Standard EN455 (Parts 1, 2 & 3) & EN374-1, 2 & 3:2003. Manufactured in a facility holding the ISO 9001:2008 and 89/686/EEC Article 11b certifications. Processed in an NEBB certified ISO Class 5 cleanroom environment.

Packaging

Double PE bagged

Sterility

Non-Sterile

Shelf life

5 years from date of manufacture

Cleanroom class

ISO Class 4+

Handling & storage

Store in a dry cool place (<40°C) away from direct sunlight



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

sales@cleanroomshop.com
Fax: + 44 (0) 1524 813022

Disposal

Recommended for single use



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

sales@cleanroomshop.com
Fax: + 44 (0) 1524 813022