

Helapet SteriClean IPA Alcohol Spray Sterile 0.9l - 1 case



Why Choose Helapet SteriClean IPA Alcohol Spray Sterile 0.9l - 1 case

Sold in cases of 6

Sterile, WFI, ready to use, 0.9L

70% v/v isopropyl alcohol (IPA) blended with 30% WFI

Double bagged & 0.2 micron filtered

Ergonomic and lightweight bottle

Bag in bottle closed system ensures sterility of contents

Helapet code: XX6800

A dangerous goods shipping surcharge will be applied to this product for shipments leaving the UK mainland

More Info

Sterile IPA Spray

SteriClean® Alcohol Spray uses an optimal mix of 70% v/v isopropyl alcohol (IPA) blended with 30% purified water for injection (WFI) and 0.2 micron filtered for use in aseptic areas.

Contained in an ergonomic and lightweight bottle design, Helapet sterile trigger sprays use a proven closed bag in bottle spray assembly to maintain sterility of contents by restricting the ingress of bacteria.

Format: 0.9L

Case of 6

Sterile, double bagged

Shelf life: 1 year from date of manufacture

Bag in bottle closed system ensures sterility of contents

Long length, non-venting trigger



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

sales@cleanroomshop.com
Fax: + 44 (0) 1524 813022

Type

Spray Gun (Trigger Spray)

Lead time (if out of stock)

5-7 working days

Pack Size

6

Materials

70% v/v isopropyl alcohol (IPA) blended with 30% purified water for injection (WFI)

Size

0.9L

Water grade

Water for injection (WFI)

Applications

Cleanrooms

Technical data

Batch control

Each SteriClean® Alcohol Spray is double bagged in an ISO Class 7 cleanroom and gamma sterilised to no less than 25 kGy, suitable for efficient aseptic transfer into cleanrooms and isolator environments.

Packaging

Double bagged

Sterility

Sterile

Shelf life

2 years from date of manufacture

Dilution

Ready to use

MSDS

Available on request

Cleanroom class

Cleanroom

Dangerous Goods

Yes



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

sales@cleanroomshop.com
Fax: + 44 (0) 1524 813022



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

sales@cleanroomshop.com
Fax: + 44 (0) 1524 813022