

Klerwipe 70/30 Denatured Ethanol Sterile Tub Wipe (Case 15)



Why Choose Klerwipe 70/30 Denatured Ethanol Sterile Tub Wipe (Case 15)

Sold in cases of 15

For bulk quantities that exceed the above pricing structure phone +44(0) 1524 813022

Binderfree, hydroentangled polyester/cellulose

Wipe size: 200x230mm

Pack size: 200 wipes

Wipe grade: 22gsm

70% IMS Tub Wipe

Ecolab code: 3079200

A dangerous goods shipping surcharge will be applied to this product for shipments leaving the UK mainland

More Info

Sterile low particulate wipes impregnated with 70|30 Sterile Denatured Ethanol.

These wipes are binderfree, hydroentangled polyester/cellulose

Wipe size: 200x230mm

Pack size: 200 wipes

Wipe grade: 22gsm

70% IMS Tub Wipe

Suggested Usage:

Grade A: Very Suitable

Grade B: Very Suitable

Grade C/D: Suitable

Unclassified: Suitable

The Klertear packaging is easy to open even when wearing gloves

Ecolab: 3079200

Please note: A dangerous goods surcharge may apply - more information can be found on our [delivery information](#) page.

Product Icon Attributes

Suitable for GMP A/B/C/D, Double bagged, Sterile, 70/30 denatured ethanol

Brand

Ecolab

Type

Tub wipes

Lead time (if out of stock)

10-15 days

Pack Size

Case of 15

Construction/Material

Non Woven

Size

200 X 230mm

Applications

Suitable for use on all types of hard surfaces in cleanrooms and controlled environments, isolators, LAF's, equipment exteriors and production vessels.

Country of origin

United Kingdom

Technical data

Available in two options:
22gsm (KW35) - 200 wipes
42gsm (KW37) - 100 wipes
Wipe size: 200x230mm

Sterility

Sterile

Shelf life

2 years

Efficacy

Efficacy
Bactericidal: EN 1276
Fungicidal: EN 1650
Bactericidal and Fungicidal: EN 13697

Handling & storage

Store in cool areas

Dangerous Goods

Yes



UKAS is a member of Registrar of Standards (Holdings) Ltd.



UKAS is a member of Registrar of Standards (Holdings) Ltd.



UKAS is a member of Registrar of Standards (Holdings) Ltd.

www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589