



## Why Choose Klerwipe Sterile Dry Wipe 10 X 10 packs

Sterile low lint, low particulate dry wipe

300mm x 300mm

10 packs x 10 wipes

Double bagged

Packed in a Grade C cleanroom

Terminally sterilised by gamma irradiation

Alternative to CCD1 - Clear Clean Dry-Wipes Sterile

Ecolab code: 3035510

## More Info

Sterile low particulate dry wipe

300mm x 300mm

10 packs x 10 wipes

Double bagged

Klertear - easy to open even when wearing gloves

Packed in a Grade C cleanroom

Terminally sterilised by gamma irradiation

68 gsm binder free hydroentangled cellulose / polyester blend wipes

Sterilised by gamma irradiation

Ecolab: 3035510

## Brand

Ecolab

## Type

Dry Wipe

