

# DuPont™ Tyvek® Isoclean® Sterile Hooded Coverall - Case of 25



## Why Choose DuPont™ Tyvek® Isoclean® Sterile Hooded Coverall - Case of 25

This Tyvek® Isoclean® coverall is a clean, processed, sterile garment, suitable for controlled environments.

### Limited stock

Sold in cases of 25

Suitable for use in GMP class A/B (ISO Class 5) clean rooms

Clean-processed and sterilised

Certificates of sterility available

Sterile packed, double bagged

Case size: 25

Size: M, XL, XXL remaining

DuPont™ code: IC0105SWHCS

## More Info

### EXTRA PROTECTION THANKS TO THE INNOVATIVE DESIGN

DuPont™ Tyvek® IsoClean® hooded coverall with attached slip-retardant shoe covers, model IC 105 S WH CS. Available in sizes MD to 3X. Clean-processed and gamma-sterilized. Serged internal seams. Thumb loops. Elastication at wrists, ankles and face. Front zipper closure.

### GOOD BALANCE BETWEEN PROTECTION, DURABILITY AND COMFORT

Tyvek® IsoClean® delivers an ideal balance of protection, durability and comfort. Made of high density polyethylene using a patented flash spinning process. Tyvek® IsoClean® provides an inherent barrier to particles, microorganisms and non-hazardous water-based light liquid splash. Tyvek® IsoClean® is also breathable and exceptionally low linting.

Tyvek® IsoClean® (option codes CS, DS and MS) garments and accessories are clean-processed to maximize cleanliness and have been sterilized by gamma-irradiation. Garments are folded to aid aseptic donning and individually packed in an



[www.cleanroomshop.com](http://www.cleanroomshop.com)

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)

Tel: + 44 (0) 1524 813022

Fax: + 44 (0) 1524 811589

XL

XXL

sterile Tyvek® IsoClean® are typically used in cleanrooms within the biotech, pharmaceutical, medical device manufacturing, food processing, cosmetics industries as well as in other critical or controlled environments.

Clean-processed and sterilised by gamma-irradiation to SAL of  $10^{-6}$  (ISO 11137-1)

Full traceability on all sterilized apparel with [certificates of sterility](#) available

Suitable for use in GMP class A/B (ISO Class 5) clean rooms

This garment is not CE certified as chemical protective clothing

## PRODUCT ORDERING

Size: (THCS1MB) Medium, (THCS1XLB) XL, (THCS1XXLB) XXL

Colour: White

Case quantity: 25

DuPont™ code: IC0105SWHCS

## Product Icon Attributes

Suitable for GMP A/B/C/D, Individually wrapped, Sterile, M-3XL

## Brand

DuPont™ Tyvek®

## Type

Hooded coverall

## Lead time (if out of stock)

3-5 days

## Pack Size

25

## Materials

Spun bonded polyethylene, non woven material

## Construction/Material

Tyvek®

## Colour

White

## Applications

Controlled environments, up to ISO Class 4 / Grade A/B aseptic cleanroom. typically used in cleanrooms within the biotech, pharmaceutical, medical device manufacturing, food processing, cosmetics industries as well as in other critical or

## Technical data

Size Chart

Body measurements in cm (inches)

Size, Chest Girth, Body Height

M: 92 100 135" 20" 1.68 1.76 1.65" 5'0"



[www.cleanroomshop.com](http://www.cleanroomshop.com)

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)

Tel: + 44 (0) 1524 813022

Fax: + 44 (0) 1524 811589

## Sterility

Sterile

## Shelf life

5 year shelf life from date of manufacture

## CE Certified

Yes

## Handling & storage

Store in a dry location and avoid direct contact with sunlight for prolonged periods.

## Disposal

Tyvek® does not contain halogens. When incinerated at high temperature Tyvek® yields only water and carbon dioxide.

## Warning

Stay away from heat and flame

## Antistatic properties

Electrostatic protection according to EN 1149-1.

## Washing

Washing and laundering is not recommended as it impacts on the performance of the materials.



[www.cleanroomshop.com](http://www.cleanroomshop.com)  
Tel: + 44 (0) 1524 813022

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)  
Fax: + 44 (0) 1524 811589