



Why Choose Klerwipe Neutral Detergent Pouch Wipe

Unit size: 300 (20 packs x 15 wipes)

Wipe size: 20 x 20cm

Impregnated with Klerclean Neutral Detergent

For oil and residue removal

Easy to open, even when wearing gloves

Price shown is per unit

Formerly - Klerwipe-CR Sterile Neutral Detergent Low Particulate Pouch Wipe

Ecolab code: 3079250

More Info

Unit size: 300 (20 packs x 15 wipes)

Wipe size: 20 x 20cm

Impregnated with Klerclean Neutral Detergent in a convenient resealable pouch

The wipes are 68gsm, binder free, hydroentangled polyester/cellulose blend wipes

Suitable for rotational cleaning

Double bagged with Klartear feature

Ecolab: 3079250

Product Icon Attributes

Sterile

Brand

Ecolab



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589

Type

Pack/Pouch wipe

Lead time (if out of stock)

4 weeks

Pack Size

300 wipes

Construction/Material

Non Woven

Size

200 x 200mm

Applications

Suitable for use on all types of hard surfaces in cleanrooms and controlled environments, isolators, LAF's, equipment exteriors and production vessels.

Country of origin

United Kingdom

Sterility

Non-Sterile

PPE Requirements

Ensure that there is sufficient ventilation to the area. Nitrile, PVC or Neoprene gloves for hand protection. Eye glasses for eye protection and protective clothing for overall protection.

Handling & storage

Store in cool well ventilated areas. Keep container tightly closed. Must only be kept in original packaging.



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589