



Why Choose Azo Universal Spray 500ml CE - Case of 12

Case of 12

Active technology, "one-step" cleaning & disinfection spray

For non-invasive medical devices

Vernacare (formerly Synergy Health) code: AM500SCE

Formerly known as Azo Max

Meets Virucidal standard EN 14476

More Info

One-step cleaning & disinfection for multi surface and medical devices

Based on a highly effective formulation of quaternary ammonium compounds (QACs) and polymeric biguanide (PHMB) working together with other ingredients to deliver a potent and fast acting disinfectant - non-irritant and non-explosive.

Active technology, "one-step" cleaning & disinfection spray

For non-invasive medical devices

CE mark - bottle

Effective bacterial action against MRSA, E.coli and a diverse range of other micro-organisms

Ideal for use in cleanrooms and critical environments

99.999% effectiveness against pathogens, validated by European Standards e.g. EN 1276, EN13727, EN 14476, EN 1040, EN1650.

One stop cleaning and disinfection - saves you time and money



www.cleanroomshop.com

info@cleanroomshop.com

Tel: + 44 (0) 1524 813022

Fax: + 44 (0) 1524 811589

Broad spectrum antimicrobial efficacy - effective against a wide range of micro-organisms.

Chlorine, glutaraldehyde and alcohol free - non-tainting and non-corrosive.

Kind to surfaces and skin friendly - gentle on surfaces and kind to your skin.

Formerly known as Azo Max

Meets Virucidal standard EN 14476

Meets Bactericidal standard EN 1276

Meets Fungicidal standard EN 16615

Product Icon Attributes

For ISO cleanrooms, Non-sterile, Broad spectrum activity, Cleans and disinfects

Brand

Vernacare

Type

Spray Gun (Trigger Spray)

Lead time (if out of stock)

10-15 days

Pack Size

Case of 12

Materials

Quaternary ammonium compounds (QACs) and polymeric biguanide (PHMB)

Size

500ml

Applications

Designed for use in: Hospitals, Nursing and Care Homes, Laboratories, Dental and GP Surgeries, Nurseries and schools

Areas for use: Local work surfaces, Patient treatment areas, Bed frames and mattresses, Examination couches, Keyboards, Phones, Medical devices class II A

Technical data

99.999% effectiveness against pathogens, validated by European Standards e.g. EN 1276, EN13727, EN 14476, EN 1040, EN1650.

Meets Virucidal standard EN 14476

Meets Bactericidal standard EN 1276

Packaging

Not bagged



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589

Long lasting residual efficacy - antimicrobial properties present for up to 24 hours after application.
Broad spectrum antimicrobial efficacy - effective against a wide range of micro-organisms.
Chlorine, glutaraldehyde and alcohol free - non-tainting and non-corrosive.
Kind to surfaces and skin friendly - gentle on surfaces and kind to your skin.

Sterility

Non-Sterile

Shelf life

2 year shelf life from the point of manufacture

Dilution

Ready to use

MSDS

Full MSDS available on request.

CE Certified

1

PPE Requirements

Non-flammable ☐ no specialist PPE required

Handling & storage

Non-flammable - no specialist storage required



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589