

Klerwipe Neutral Detergent Mop Wipe - 10 Pouches



Why Choose Klerwipe Neutral Detergent Mop Wipe - 10 Pouches

Unit size: 200 (10 x 20 wipe pouches)

Wipe size: 42 x 25cm

Impregnated with neutral detergent & sterile

Compatible with Product Codes: MS3, MS4

Price shown is per unit

Ecolab code: 3079260

More Info

Advantages

Effective - Sterile low particulate mop wipes impregnated with neutral detergent. Klerclean-CR is particularly effective on oily and sticky spillages or following general maintenance.

It is a non-ionic sterile detergent, a mix of deionised water and a fatty alcohol ethoxylate.

The detergent is 0.2 micron filtered then filled and bagged in a Grade C cleanroom before being terminally sterilised by gamma irradiation. A sterile ready-to-use detergent mop wipe, in a peelable pack, for large surfaces.

Pouch - Each pouch contains 20 wipes to reduce wastage during a manufacturing session. The double bagging features Klertear which makes it easy to open even when wearing gloves. The pouch is easily resealable to prevent the wipes drying out and does not itself generate any particles. Wipes are available as 68 gsm binder-free, hydroentangled cellulose / polyester.

Product ordering

Format: Impregnated Pouch Wipes 420 x 250 mm

Case quantity: 10 x 20 wipes

Ecolab code: 3079260

[Product datasheet \(PDF\)](#)

Product Icon Attributes

Sterile, Neutral detergent, Polyester / cellulose, 10 packs of 20 wipes

Brand

Ecolab

Type

Mop Wipe

Lead time (if out of stock)

4 weeks

Pack Size

200 (10x20)

Construction/Material

Non Woven

Size

42cmx25cm

Colour

White

Country of origin

United Kingdom

Sterility

Sterile



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589