



Why Choose Klercide 70/30 IPA Filtered Capped 4 x 5L

Filtered, DI, ready to use, 5 litre

Non-sterile

Cleanroom filtered 70% IPA

70% alcohol blended with deionised water

Double bagged & 0.2 micron filtered

Great value 4 x 5L capped bottle

Suitable for a high number of applications in Grades C/D/Unclass areas

Ecolab code: 3078420

A dangerous goods shipping surcharge will be applied to this product for shipments leaving the UK mainland

More Info

Filtered Klercide 70|30 IPA is a blend of 70% v/v Isopropyl Alcohol blended with deionized water. The alcohol blend is 0.2 micron filtered, filled and double bagged in a Grade C (ISO Class 7) cleanroom.

Advantages

Easy to use - Convenient range of formats available; a comfortable, easy to use 1 litre trigger spray and a large 5 litre capped version.

More environmentally friendly - The product contains no propellants so does not require any special disposal procedures.

Product ordering

Format: 5 Litre Capped

Case quantity: 4 x 5L

Ecolab code: 3078420

[Safety datasheet \(PDF\)](#)

Brand

Ecolab

Type

Capped Bottle

Lead time (if out of stock)

10-15 days

Pack Size

1

Materials

70% isopropyl alcohol blended with deionised water

Size

5L

Colour

Colourless

Odour

Alcohol-like

Clarity

Clear

Water grade

DI

Applications

Ready to use so can be sprayed directly onto the surface to be sanitised. Always ensure the whole surface is in contact with the liquid. For optimum results wipe the surface to physically remove any contaminants.

Production

Klercide 70 | 30 IPA blended with deionized water is 0.2 micron filtered then filled and double bagged in a Grade C cleanroom

Country of origin

United Kingdom

Technical data

Suitable for a high number of applications in Grades C/D/Unclass areas - please enquire for full details

Packaging

Double Bagged

Sterility

Non-Sterile

Shelf life

2 years from date of manufacture

Efficacy

CEN Test Data available on request. Certificates of analysis,

irradiation, quantitative endotoxin levels < 0.25EU/ml and sterility available

Dilution

Ready to use

MSDS

Available to download on the Cleanroomshop product page

PPE Requirements

Ensure that there is sufficient ventilation to the area. Nitrile, PVC or Neoprene gloves for hand protection. Eye glasses for eye protection and protective clothing for overall protection.

Handling & storage

Store in cool well ventilated areas. Keep container tightly closed. Must only be kept in original packaging.

Dangerous Goods

Yes

Instructions for use

Ready to use so can be applied directly onto the surface to be sanitised. Always ensure the whole surface is in contact with the liquid. For optimum results wipe the surface to physically remove any contaminants.



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589