



Why Choose KIMTECH A7 P+ Lab coat - Case of 15

For protection against a wide variety of light chemical splash events, KIMTECH SCIENCE* A7 P+ white lab coats have your processes covered.

Cat III, Type 6 (PB)

Limited chemical splash protection, limited chemical splash protection

High levels of particle and bacterial protection (with filtration efficiency of 99%)

Cloth-like low lint fabric

Three layers of fabric offering strength, durability, and protection

Double-bagged for ISO 7, Grade C areas

Case size: 15

Size: S - 2XL

LIMITED STOCK!

More Info

KIMTECH SCIENCE* A7 P+ white lab coats feature elasticated cuffs, thumb-loops, extra length arms and high mandarin collars to enhance all-round coverage and security.

Kimtech™ A7 P+ Laboratory Coats are disposable laboratory wear designed to maximise safety, comfort and compliance. The multi-layered coats offer secure coverage against light chemical splash events and other contaminants, for seamless protection when and where it counts. The latex-free and non-sterile coats are suitable for a variety of applications in environments that demand basic contamination control. The multi-layered design has high levels of particle and bacterial protection (with filtration efficiency of 99%) and is formed of non-woven SMS fabric with triple stitched seams and abrasion-resistant, film coated PP. The material has a non-slip outer layer and a cloth-like feel, while also being very breathable, reducing the risk of heat stress and making the lab coats suitable for daily use. Kimtech™ A7 P+ white lab coats also feature elasticated cuffs, thumb-loops, extra length arms and high mandarin collars to enhance all-round coverage and security. Both the chest and wrist areas are fully covered, and lab gloves can be secured over the coat, giving excellent coverage. The secure construction makes the coats suitable for demanding processes in pharmaceutical, biochemistry, biomedical, cosmetics, forensics, and other

Available Options

Small

Medium

Large

XL

XXL

environments. Silicon-free and low-lint, Kimtech™ A7 P+ disposable lab coats are individually-packed and double-bagged, making them suitable for ISO 7 (Grade C) areas. The laboratory coats are available in sizes S to 2XL.

Assured Compliance:- PPE Cat III according to (EU) Regulation 2016/425 :Type 6 (PB) limited chemical splash protection- EN1149-1 Antistatic clothing - Tensile strength EN ISO 13934-1 Class 1 of 3 - Surface resistivity EN 1149-1 < 5 x 10¹⁰ ohm- Particle shedding test (Helmke Drum) IEST-RP-CC003.3 Category II- Bacterial Filtration Efficiency (3.0 µm) ASTM F2100 >99%

High quality disposable laboratory coats with triple stitched seams and breathable, multi-layered material.

Enhanced all-round protective features are included, such as high mandarin collars, thumb-loops and extra-length arms to fit lab gloves over coats.

The breathable coat material includes high levels of particle and bacterial protection (with filtration efficiency of 99%).

The coats are individually packed making them suitable for ISO 7 (Grade C) areas.

Online ordering

Manufacturers code: A7 P+

Colour: White

Pack size: 15

Size: S (KCLC2S), M (KCLC2M), L (KCLC2L), XL (KCLC2XL), XXL (KCLC2XXL)

LIMITED STOCK!

Product Icon Attributes

Suitable for ISO 7-8, Double bagged, Non-sterile, S-2XL, Cat III Type PB (6B)

Brand

Kimberly Clark

Type

Lab coat

Lead time (if out of stock)

10-15 Days

Materials

Non-woven SMS Fabric

Applications

The coats are individually packed making them suitable for ISO 7 (Grade C) areas.

Technical data

PPE Cat III according to (EU) Regulation 2016/425 :Type 6 (PB)
limited chemical splash protection- EN1149-1 Antistatic
clothing - Tensile strength EN ISO 13934-1 Class 1 of 3 -
Surface resistivity EN 1149-1 99%

Sterility

Non-Sterile



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589