

# Klercide Neutral Detergent Sterile Spray 1L - Case of 6



## Why Choose Klercide Neutral Detergent Sterile Spray 1L - Case of 6

Case of 6

Sterile, ready to use, 1 litre

Removes contamination

Double bagged & 0.2 micron filtered

Suitable for rotational cleaning

Very suitable for a number of applications in Grades A/B/C/D areas

Ecolab code: 3078790

## More Info

Klercide Neutral Detergent is a non-ionic, sterile neutral detergent with a pH between 5.5–7.5. It is a mix of deionized water and a fatty alcohol ethoxylate. It is 0.2 micron filtered, filled and double bagged in a Grade C (ISO Class 7) cleanroom before gamma irradiation at no less than 25 kGy.

It is designed to remove contamination in critical areas particularly cleanrooms, isolators, production equipment and all nonporous surfaces. It is particularly effective on oily and sticky spillages or following general maintenance especially where lubricants have been applied.

The trigger spray incorporates the patented SteriShield Delivery System and an adjustable trigger. The unique way in which the SteriShield Delivery System works, as a closed system, ensures that the sterility of the contents is preserved throughout use. This is especially important for detergents which have no inherent biocidal activity and the potential consequences, should microbiological contamination occur, are particularly serious.

The ready to use product is available as a 1 litre trigger spray and a 5 litre capped product.

## Advantages

Very effective - Klercide Neutral Detergent is particularly effective on oily and sticky spillages or following general maintenance. It is suitable for surface cleaning in gassing isolators.

Easy to use - It is classified as non-corrosive so it can be used in most cleanroom applications. The fully adjustable trigger spray allows liquid to be dispensed as either a jet or a spray, which is especially useful in critical areas.

Suitable for direct product contact - Klercide Neutral Detergent is a sterile, non-ionic, low foaming neutral detergent, with low residues, which are easily removed by rinsing.

## Product ordering

Format: 1 Litre SDS Trigger Spray

Case quantity: 6 x 1L

Ecolab code: 3078790

[Safety datasheet \(PDF\)](#)

## Product Icon Attributes

Suitable for GMP A/B/C/D, Double bagged, Sterile, Neutral detergent

## Brand

Ecolab

## Type

SteriShield delivery system

## Lead time (if out of stock)

10-15 days

## Pack Size

6

## Materials

A blend of deionised water and a fatty alcohol ethoxylate

## Size

1L

## Colour

Colourless

## Odour

Faint

## Clarity

Clear

## Applications

Very suitable for a number of applications in Grades A/B/C/D areas - please enquire for full details

## Country of origin

United Kingdom

## Sterility

Sterile

## Dilution

Ready to use

## Cleanroom class

GMP Grade A, GMP Grade B, GMP Grade C, GMP Grade D

## Handling & storage

Store in cool well ventilated areas. Keep container tightly closed.

Must only be kept in original packaging.

## Production

Klercide Neutral Detergent is filtered to 0.2 micron and filled in a Grade C (ISO Class 7) cleanroom.

## Packaging

Double Bagged

## Shelf life

2 years from date of manufacture

## MSDS

Available to download from the Cleanroomshop product page

## PPE Requirements

Ensure there is sufficient ventilation of the area.

Hand protection: Protective gloves.

Eye protection: Safety glasses.

Skin protection: Protective clothing with elasticated cuffs and closed neck.

## Delivery system

The way in which the SteriShield Delivery System works, as a closed system, ensures that the sterility of the product is preserved throughout use.



[www.cleanroomshop.com](http://www.cleanroomshop.com)  
Tel: + 44 (0) 1524 813022

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)  
Fax: + 44 (0) 1524 811589