



## Why Choose PurQuip® Uno Cleanroom Mopping System

Compact stainless steel cleaning trolley system with wringer

2 stainless steel buckets with capacities of 8.5 litres at the front

Ideal for ISO Class 4-8 & GMP A-B high end pharma areas

Uno 200: no rear buckets (2121200)

Uno 211: with rear buckets (2121203)

## Available Options

Uno 200

Uno 211

## More Info

PurQuip® Uno | Stainless steel cleaning trolley

Top quality trolley system for cleanroom cleaning

Completely made of stainless steel

Double bucket system (2 x 8.5 litre) with optional rear buckets

Mirror polished premium surface

Ideal for high end pharma production areas

Durable, autoclavable, disinfectant resistant

Very compact and manoeuvrable (small space requirement)

Sturdy due to optimal weight distribution

Conforms to GMP, DIN ISO and VDI

Castors have grease-free bearings and are trackless

Compatible with PurMop® mops & PurQuip® accessories

Recommended for up to cleanroom class ISO 4

## Options available:

MSHPQ14 - Uno 200 mop trolley with press - 2 x 8.5L bucket

system without extra buckets (2121200)

MSHPQ15 - Uno 211 mop trolley with press - 2 x 8.5L bucket system with 2 x 8.5L buckets at rear (2121203)

## Product Icon Attributes

Suitable for ISO 4-8, Auto clavable, Trolley system, Stainless steel

## Brand

Hydroflex

## Lead time (if out of stock)

Up to 6 weeks

## Pack Size

1

## Materials

Stainless steel

## Construction/Material

Stainless steel

## Applications

Recommended for up to cleanroom class ISO 4 [gMP A/B]

## Country of origin

Germany

## Technical data

Type Uno 100 | Single bucket system with flat-press - MSHPQ1

1 stainless steel flat-press type

1 stainless steel container, type (26 Ltr.) on the front

Size: 46 x 58 x 89cm

Type Uno 111 | Single bucket system with flat-press and extra containers - MSHPQ13

1 stainless steel flat-press

1 stainless steel container, (26 Ltr.) on the front

2 stainless steel containers, (8.5 Ltr.) at the rear

2 stainless steel support brackets (can be hung on the rear)

Size: 56 x 58 x 89cm

Type Uno 200 | Double bucket system with flat-press - MSHPQ14

1 stainless steel flat-press

2 stainless steel containers, (8.5 Ltr.) on the front

Size: 46 x 58 x 89cm

Type Uno 211 | Double bucket system with flat press & additional container - MSHPQ15

1 stainless steel flat-press

2 stainless steel containers, (8.5 Ltr.) on the front

2 stainless steel containers, (8.5 Ltr.) at the rear

2 stainless steel support brackets (can be hung on the rear)

Size: 56 x 58 x 89cm

## Sterility

Non-Sterile

## Castors

Conforms to GMP, DIN ISO and VDI



[www.cleanroomshop.com](http://www.cleanroomshop.com)  
Tel: + 44 (0) 1524 813022

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)  
Fax: + 44 (0) 1524 811589