


SHIELDskin XTREME™ Sterile White Nitrile 330 DI+ Gloves



Why Choose SHIELDskin XTREME™ Sterile White Nitrile 330 DI+ Gloves

Powder-free triple DI washed hand-specific sterile 33cm nitrile gloves with textured palm and fingers.

Nitrile

Cleanroom packed, sterile, double bagged

Length: 330mm

Sizes: 5.5-10

Pack size: 1 pair/pouch, 20 pouches/bag, 10 bags/case (200 pairs)

Colour: white

PPE Cat III, AQL 0.65

Virus resistant according to ISO 16604 Procedure B & ISO 374-5:2016

Suitable for use in EU GMP Grade A/B (ISO Class 5) cleanrooms

SHIELD Scientific code: 69 876X Series

More Info

Powder-free triple DI washed hand-specific sterile 33cm nitrile gloves with textured palm and fingers.

Nitrile

Cleanroom packed, sterile, double bagged

Length: 330mm

Sizes: 5.5-10

Pack size: 1 pair/pouch, 20 pouches/bag, 10 bags/case (200 pairs)

Colour: white

PPE Cat III, AQL 0.65

Virus resistant according to ISO 16604 Procedure B & ISO 374-5:2016

Suitable for use in EU GMP Grade A/B (ISO Class 5) cleanrooms

SHIELD Scientific code: 69 876X Series

Category III PPE glove - PPE Regulation (EU) 2016/425

Extensively tested for chemical permeation according to EN 16523-1:2015

Available Options

6.5
7.0
7.5
8.0
8.5
9.0

Static dissipative according to EN 1149-1/2/3 & 5

Low endotoxin content at <20 EU/pair (EN 455-3:2015)

Longer than normal length to provide enhanced grip pf garment sleeve

Ultra clean with specification limit of <1200 particles per cm² @≥0.5µm

Tested to ASTM D6978-05 to demonstrate suitability for handling cytotoxic drugs

Lot-specific certificate of conformance detailing barrier performance, extractables, particles and endotoxin downloadable from website

Highest level of micro-organism resistance (Level 3, EN 374-2:2014)

Technical Data Pack available on request

Chemical risk (ISO 374-1:2016 + A1:2018 Type B JKP)

Product Icon Attributes

Nitrile, Sterile, Cat III, Double Bagged

Brand

SHIELD Scientific

Lead time (if out of stock)

Enquire for lead time

Pack Size

1 pair/pouch, 20 pouches/bag, 10 bags/case

Colour

White

Country of origin

Malaysia

Sterility

Sterile

Shelf life

5 years



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589