



Why Choose SHIELDskin XTREME™ Sterile Latex 400 DI+

Natural rubber latex and fully textured gloves

Sterile, double bagged

Powder-free, triple DI washed, hand-specific

Length: 400mm

Sizes: 5.5-10.0

Pack size: 1 pair/pouch, 20 pouches/bag, 8 bags/case (160 pairs)

Colour: Natural

PPE Cat III, AQL 0.65

SHIELD Scientific code: 69 577X Series

Virus resistant according to ISO 16604 Procedure B & **ISO 374-5:2016**

Available Options

More Info

69 577X series

STERILE: Biocontamination control

Natural rubber latex and fully textured gloves

Sterile cleanroom glove

Category III PPE glove (PPE Regulation 2016/425)

Complex Design - For mortal & irreversible risks

Powder-free natural latex glove

Hand-specific

400 mm / 0.18 mm (EN 420:2003+A1:2009)

Biological risk (ISO 374-5:2016 VIRUS)

AQL 0.65 (EN 374-2:2014 Level 2)

Chemical risk (ISO ISO 374-1:2016+A1:2018 - Type B KPT)

Waterproof and for low chemical protection

Tested for chemical permeation (EN 16523-1:2015+A1:2018)

Typical particle levels: less than 1100 per cm² ≥ 0.5µm

Type I hypersensitivity & type IV hypersensitivity reduced

5.5	Suitable for use in EU GMP Grade A/B (ISO Class 5) cleanrooms
6.0	
6.5	Lot-specific certificate of conformance detailing barrier performance, extractables, particles and endotoxin downloadable from website
7.0	
7.5	Technical Data Pack available on request
8.0	Virus resistant according to ISO 16604 Procedure B & ISO 374-5:2016
8.5	Options:
9.0	SHIELDskin Xtreme Sterile Latex 400 DI+ 5.5 - SSXSL2155
10.0	SHIELDskin Xtreme Sterile Latex 400 DI+ 6.0 - SSXSL2160
	SHIELDskin Xtreme Sterile Latex 400 DI+ 6.5 - SSXSL2165
	SHIELDskin Xtreme Sterile Latex 400 DI+ 7.0 - SSXSL2170
	SHIELDskin Xtreme Sterile Latex 400 DI+ 7.5 - SSXSL2175
	SHIELDskin Xtreme Sterile Latex 400 DI+ 8.0 - SSXSL2180
	SHIELDskin Xtreme Sterile Latex 400 DI+ 8.5 - SSXSL2185
	SHIELDskin Xtreme Sterile Latex 400 DI+ 9.0 - SSXSL2190
	SHIELDskin Xtreme Sterile Latex 400 DI+ 10.0 - SSXSL2110

Brand

SHIELD Scientific

Lead time (if out of stock)

.

Pack Size

1 pair/pouch, 20 pouches/bag, 8 bags/case (160 pairs)

Construction/Material

Latex

Country of origin

Malaysia

Sterility

Sterile

Shelf life

5 years



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589