



Why Choose SHIELDskin XTREME™ Sterile White Nitrile 400 DI+ Gloves

Nitrile gloves with textured palm and fingers

Sterile

Powder-free, triple DI washed, hand-specific

Length: 400mm

Sizes: 5.5-10.0

Pack size: 1 pair/pouch, 20 pouches/bag, 8 bags/case (160 pairs)

Colour: white

PPE Cat III, AQL 0.65

SHIELD Scientific code: 69 877X series

Virus resistant according to ISO 16604 Procedure B & **ISO 374-5:2016**

More Info

69 877X series

STERILE: Biocontamination control

Sterile cleanroom glove

Category III PPE glove (PPE Regulation (EU) 2016/425)

Complex Design - For mortal & irreversible risks

Powder-free white nitrile glove

Hand-specific

400 mm / 0.15 mm (EN 420:2003+A1:2009)

AQL 0.65 (EN 374-2:2014 Level 3)

Biological risk (ISO 374-5:2016 VIRUS)

Chemical risk (ISO 374-1:2016 + A1:2018 Type B KPT)

Tested for chemical permeation (EN 16523-1:2015+A1:2018)

Typical particle levels: less than 1000 per cm² ≥ 0.5µm

Tested according to EN 1149-1-2-3 & 5

Type I hypersensitivity eliminated - Type IV hypersensitivity

Available Options

5.5
6.0
6.5
7.0
7.5
8.0
8.5
9.0
10.0

reduced

Suitable for use in EU GMP Grade A/B (ISO Class 5) cleanrooms

Lot-specific certificate of conformance detailing barrier performance, extractables, particles and endotoxin downloadable from website

Technical Data Pack available on request

Virus resistant according to ISO 16604 Procedure B & **ISO 374-5:2016**

Product Icon Attributes

Nitrile, Sterile, Category III, Virus Resistant

Brand

SHIELD Scientific

Lead time (if out of stock)

Enquire for lead time

Pack Size

1 pair/pouch, 20 pouches/bag, 8 bags/case (160 pairs)

Construction/Material

Nitrile

Colour

White

Country of origin

Malaysia

Sterility

Sterile

Shelf life

5 years



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589