

Agma Sterile Zyceine in WFI wipes 68gsm 10 x 5 wipes



Why Choose Agma Sterile Zyceine in WFI wipes 68gsm 10 x 5 wipes

This product has a temporary change in packaging - please see below

Optimised fluid dosing in conjunction with the highly absorbent 68 gsm polyester/ cellulose blend enables the wipe to hold sufficient fluid for in-use requirements without drying out too quickly

To be used in all grades of cleanroom (ABCD)

Sterile packed in a zipock sealed plastic bag

5 Wipes per pouch (Box dimensions 29x22x17.5cm)

10 pouches containing 5 wipes per pouch (50 wipes per case)

The goods are classed as non-haz

More Info

Optimised fluid dosing in conjunction with the highly absorbent 68 gsm polyester/ cellulose blend enables the wipe to hold sufficient fluid for in-use requirements without drying out too quickly

Sterile packed in a zipock sealed plastic bag.

Supplied quarter folded for easy wiping convenience

Each wipe is 28 x 38cm, allowing a large surface area to be covered

Packed in a triple laminate 0.2 micron filtered and 25kGy gamma irradiated for sterility

Manufactured using WFI quality water guaranteeing a low endotoxin level and suitability for product contact areas

Batch coded with Detex dot to clearly indicate successful irradiation for peace of mind

2 year unopened shelf life"

Hypochlorous acid solution (Aluminium foil packaging)

Temporary change of packaging from 21st November

The manufacturer has temporarily changed the packaging for this product as its primary source of foil is not available. The secondary source of foil is of the same specification but is white instead of silver. If you would like any more information on the change, including the results of sterilisation tests, please [contact us](#).

Product Icon Attributes

Suitable for GMP A/B/C/D, Sterile, WFI, Zyceine

Brand

AGMA

Lead time (if out of stock)

10-15 days

Pack Size

10 pouches containing 5 wipes per pouch (50 wipes per case).

Construction/Material

Hypochlorous acid solution

Size

29 x 22 x 17.5cm

Country of origin

United Kingdom

Sterility

Sterile

Shelf life

2 years from date of manufacturer



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589