

Agma Sterile 70% IPA in WFI 68 gsm 10 x 10 wipes



Why Choose Agma Sterile 70% IPA in WFI 68 gsm 10 x 10 wipes

Each wipe is 28 x 38cm, allowing a large surface area to be covered

To be used in all grades of cleanroom (ABCD)

Optimised fluid dosing in conjunction with the highly absorbent 68 gsm polyester/cellulose blend enables the product to hold sufficient fluid for in-use requirements without drying out too quickly

Supplied quarter folded for easy wiping convenience

Doesn't contain any chemical binders and so no residue is left behind in use

Alcohol blend is 0.2 micron filtered and 25kGy gamma irradiated for sterility

More Info

Manufactured in the UK by AGMA

Made with WFI which is EP/USP grade specification

10 wipes per pouch, 10 pouches per case (100 wipes per case)

28 x38 polycellulose wipe. (Box dimensions 29 x 22 x 17.5cm)

Manufactured by AGMA who work in accordance with cGMP guidelines

AGMA are ISO9001:2015 accredited

Utilised by NHS cleanrooms and global pharmaceutical companies

Sterile, packed in a triple laminate, re-sealable foil bag which is easily opened and closed whilst wearing gloves

Batch coded with Detex dot to clearly indicate successful irradiation for peace of mind

2 year unopened shelf life

Brand

AGMA

Lead time (if out of stock)

10-15 days

Pack Size

10 wipes per pouch, 10 pouches per case (100 wipes per case).

Construction/Material

Polycellulose

Size

28 x 38cm

Country of origin

United Kingdom

Sterility

Non-Sterile

Shelf life

2 years from date of manufacture



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589