

PurMop ICT2050 Isolator Cleaning Tool with fixed handle



Why Choose PurMop ICT2050 Isolator Cleaning Tool with fixed handle

PurMop® ICT2050 Isolator cleaning tool

Polished stainless steel

1 x fixed handle

Reach: 65cm

1 x mop pad frame

Sold individually

Recommended for Isolators, LAF workbenches, RABS, GMP cleanroom grade A-B or ISO 5-6

Compatible with Hydroflex mop pads

Hydroflex code: 2131015

More Info

PurMop® ICT2050 Isolator cleaning tool

Polished stainless steel

1 x fixed handle: 50cm

Reach: 65cm

1 x mop pad frame: 20 cm

Sold individually

Individually packaged

Stainless steel

Recyclable

Recommended for Isolators, LAF workbenches, RABS, GMP cleanroom grade A-B or ISO 5-6

Compatible with Hydroflex mop pads

Hydroflex code: 2131015



www.cleanroomshop.com

Tel: + 44 (0) 1524 813022

info@cleanroomshop.com

Fax: + 44 (0) 1524 811589

Brand

Hydroflex

Lead time (if out of stock)

Up to 6 weeks

Pack Size

1

Construction/Material

Stainless steel

Size

50 cm

Colour

Grey

Country of origin

Germany

Packaging

Individually Packed

Sterility

Non-Sterile

Shelf life

N/A

Warranty

None

Cleanroom class

ISO Class 5+



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589