

PurMop ICT2080 Isolator Cleaning Tool with Telescopic handle



Why Choose PurMop ICT2080 Isolator Cleaning Tool with Telescopic handle

PurMop® ICT2080 Isolator cleaning tool

Polished stainless steel

1 x telescopic handle

Reach: 62-93 cm

1 x mop pad frame

Sold individually

Recommended for Isolators, LAF workbenches, RABS, GMP cleanroom grade A-B or ISO 5-6

Compatible with Hydroflex mop pads

Hydroflex code: 2131024

More Info

PurMop® ICT2080 Isolator cleaning tool

Polished stainless steel

1 x telescopic handle: 80cm (2x40cm)

Reach: 62-93 cm

1 x mop pad frame: 20 cm

Sold individually

Individually packaged

Stainless steel

Recyclable

Recommended for Isolators, LAF workbenches, RABS, GMP cleanroom grade A-B or ISO 5-6

Compatible with Hydroflex mop pads

Hydroflex code: 2131024

Brand

Hydroflex

Lead time (if out of stock)

6-8 weeks

Pack Size

1

Construction/Material

Stainless steel

Size

40cm

Colour

Grey

Country of origin

Germany

Packaging

Individually packed

Sterility

Non-Sterile

Shelf life

n/a

Warranty

None

Cleanroom class

ISO Class 5+



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589