

PurMop ICTA2080 Isolator Cleaning Tool Aluminum Handle



Why Choose PurMop ICTA2080 Isolator Cleaning Tool Aluminum Handle

PurMop® ICTA2080 Isolator cleaning tool

Aluminium, autoclavable

1 x telescopic handle

Reach: 55-87cm

1 x mop pad frame

Sold individually

Recommended for Isolators, LAF workbenches, RABS, GMP cleanroom grade A-B or ISO 5-6

Compatible with Hydroflex mop pads

Hydroflex code: 2131028

More Info

PurMop® ICTA2080 Isolator cleaning tool

Aluminium, autoclavable

1 x telescopic handle: 80cm (2x40cm)

Reach: 55-87cm

1 x mop pad frame: 20cm

Sold individually

Individually packaged

Recyclable

Recommended for Isolators, LAF workbenches, RABS, GMP cleanroom grade A-B or ISO 5-6

Compatible with Hydroflex mop pads

Hydroflex code: 2131028

Brand

Hydroflex

Lead time (if out of stock)

Up to 6 weeks

Pack Size

5

Construction/Material

Aluminium

Size

40cm

Colour

Grey

Country of origin

Germany

Packaging

Individually packed

Sterility

Non-Sterile

Shelf life

n/a

Warranty

None

Cleanroom class

ISO Class 5+



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589