



Why Choose ecoSHIELD™ Eco Nitrile PF 250

Powder-free ambidextrous extra length (250-260 mm / 9.8"-10.2") non-sterile nitrile protective gloves.

Nitrile

Non-sterile

twinSHIELD™ double walled protection

Length: 250-260 mm

Sizes: XS-XL

Pack size: 150 gloves per dispenser - 10 dispensers per carton. (1500 gloves)

Colour: Green outer, white inner

PPE Cat III, AQL 0.25

Virus resistant according to ISO 16604 Procedure B & ISO 374-5:2016

SHIELD Scientific code: 62 512X Series

More Info

Powder-free ambidextrous extra length (250-260 mm / 9.8"-10.2") non-sterile nitrile protective gloves.

Nitrile

Non-sterile

Length: 250-260 mm

Sizes: XS-XL

Pack size: 150 gloves per dispenser - 10 dispensers per carton. (1500 gloves)

Colour: Green outer, white inner

PPE Cat III, AQL 0.25

Virus resistant according to ISO 16604 Procedure B & ISO 374-5:2016

SHIELD Scientific code: 62 512X Series

Personal Protective Equipment Category III (PPE - Complex Design) according to Regulation (EU) 2016/425.

Medical Device Class 1 (MDD) according to the Directive 93/42/EEC.

Available Options

XS

Small

Medium

Large

XL

Product Icon Attributes

Nitrile, Non-sterile, Category III, XS-XL

Lead time (if out of stock)

Enquire for lead time

Construction/Material

Nitrile

Country of origin

Malaysia

Fully compliant to the latest EU PPE norms relating to protective gloves against chemicals, micro-organisms and viruses.

twinSHIELD™ double-walled protection to offer a stronger glove and to reduce the risk of pinholes.

Chemical risk (ISO 374-1:2016 + A1:2018 Type B JKPT)

Two colours: green to make it easier to select according to the risk, combined with a soft and comfortable white interior

Brand

SHIELD Scientific

Pack Size

150 gloves per dispenser - 10 dispensers per carton

Colour

Green

Sterility

Non-Sterile



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589