



Why Choose Cleanroom Nonwoven Wipes Polyester/Cellulose 9"x9"

Non-sterile

Used for disinfection, spillages or general cleaning

45% polyester and 55% cellulose

Nonwoven - high absorbency

Triple packed & ISO 6 / Grade D Compatible

Wipe size: 230 x 230mm (9")

Pack size: 300 wipes x 10

Hydroflex code: 2183100

Product Icon Attributes

Suitable for ISO 6-8, Non-sterile, 9 x 9", 10 packs of 300 wipes

More Info

PurWipe® N1-109 Cleanroom Nonwoven Wipes Polyester / Cellulose 9" x 9"

PurWipe® N1 polyester-cellulose wipes consist of a quality, hydroentangled nonwoven material. Polyester-cellulose nonwoven

Highly absorbent material

Cost effective

Suited for general cleaning tasks

VHP-resistant, cleanroom bag with easy-tear, triple packed, for safe entry into the cleanroom

Material weight: 68g/m²

Thanks to its high absorbency and material structure, this wipe is well suited for disinfection, removal of spillages or general cleaning tasks.

Cleanroom compatible: ISO 6/7/8 GMP GRADE D

Colour: white

Case size: 3000 wipes (300 wipes / pack, 10 packs / carton)

Brand

Hydroflex

Type

Pack/Pouch wipe

Lead time (if out of stock)

Up to 6 weeks

Pack Size

3000 wipes (300 wipes / pack, 10 packs / carton)

Construction/Material

Polyester / Cellulose

Size

9" x 9"

Colour

White

Applications

Recommended for various applications such as cleaning, disinfection, liquid absorption in lower classified cleanrooms and on non-critical surfaces

Country of origin

China

Packaging

Triple packed

Sterility

Non-Sterile

Shelf life

3 years

Cleanroom class

ISO Class 6+



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589