

Sterile Cleanroom Polyester Knit Wipes 100% Polyester 12x12"



Why Choose Sterile Cleanroom Polyester Knit Wipes 100% Polyester 12x12"

Sterile

High absorbency and good general cleaning

100% Polyester

Low particle emission

ISO 4 / Grade A Compatible

Wipe size: 310 x 310mm (12")

Case size: 1000

Hydroflex code: 2183205

More Info

PurWipe® K1-212-S Cleanroom Polyester knit Wipes 100% Polyester 12"x 12"

The PurWipe® K1 is a high quality polyester cleanroom wipe with lowest particle emission which is suitable for diverse applications in critical cleanroom environments. It consists of a robust double knit which is fully decontaminated and boasts a high fluid absorption.

Robust polyester double knit. Lowest particle emission, fully decontaminated. High absorbency and good general cleaning ability. Optimal for general cleaning and disinfection on critical surfaces. VHP-resistant packing. Heat sealed edges. Material Weight 145g/m²

Cleanroom compatible: ISO 4/5/6/7/8 GMP GRADE A/B/C/D

Colour: white

Case size: 1000 wipes (10 wipes / bag, 10 bags / pack, 10 packs / carton)

Brand

Hydroflex

Type

Pack/Pouch wipe

Lead time (if out of stock)

2-3 weeks

Pack Size

1000 wipes (10 wipes / bag, 10 bags / pack, 10 packs / carton)

Construction/Material

Polyester

Size

12" x 12"

Colour

White

Applications

Recommended for cleaning and disinfection of critical surfaces in GMP and ISO cleanrooms

Country of origin


China

Sterility

Sterile

Shelf life

2 years



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589