



Sterile Polyester 2-ply Knit Wipes 100% Polyester 12"x12"



Why Choose Sterile Polyester 2-ply Knit Wipes 100% Polyester 12"x12"

Sterile

High absorption capacity

100% Polyester

Low particle emission, decontaminated

ISO 4 / Grade A Compatible

Wipe size: 310 x 310mm (12")

Case size: 800 wipes

Hydroflex code: 2183603

More Info

PurWipe® K4-312-S Sterile Polyester 2-ply knit Wipes 100% Polyester 12"x 12"

The premium cleanroom wipe PurWipe® K4 consists of a double-layered robust polyester knit with a special bonding. Thanks to its heavy construction it boasts a high absorption capacity and is suitable for routine disinfection or removal of spillages. This wipe is abrasion resistant and also applicable on rougher surfaces or stainless steel.

2-ply bonded polyester double knit. High absorption capacity. Low particle emission, decontaminated. Suitable for rough surfaces, machines, equipment. VHP-resistant packing. Ultrasonic sealed edges. Material Weight 270 g/m².

Cleanroom compatible: ISO 4/5/6/7/8 GMP GRADE A/B/C/D

Colour: white

Case size: 500 wipes (5 wipes / bag, 10 bags / pack, 10 packs / carton)

Product Icon Attributes

Suitable for GMP A/B/C/D, Sterile, 10 packs of 10x5 wipes, 12 x 12"

Brand

Hydroflex

Type

Pack/Pouch wipe

Lead time (if out of stock)

Up to 6 weeks

Pack Size

500 wipes (5 wipes / bag, 10 bags / pack, 10 packs / carton)

Construction/Material

Polyester

Size

12" x 12"

Applications

Recommended for cleaning and disinfection of surfaces, equipment and machines in higher classified cleanrooms.

Country of origin

China

Sterility

Sterile

Shelf life

2 years

Cleanroom class

ISO Class 4+, GMP Grade A



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589