

Klerwipe Sporicidal Low Residue Peroxide Mop Wipe



Why Choose Klerwipe Sporicidal Low Residue Peroxide Mop Wipe

Unit size: 200 (10 x 20 wipe pouches)

BPR compliant

Wipe size: 42 x 25cm

Impregnated with Biocide C & sterile

Can be used up to 10x without residue removal

Formerly Biocide C

Formerly - Klerwipe-CR Sterile Biocide C Low Particulate Mop Wipe

Ecolab code: 6600554

More Info

Effective - Sterile low particulate mop wipes impregnated with a very low residue blend of 6% hydrogen peroxide. The biocide blend is 0.2 micron filtered and sterile filled under Grade A laminar flow and bagged in a Grade B cleanroom. A sterile ready-to-use very low residue rotational sporicide mop wipe for large surfaces.

Pouch - Each pouch contains 20 wipes to reduce wastage during a manufacturing session. The double bagging features Klertear which makes it easy to open even when wearing gloves. The pouch is easily resealable to prevent the wipes drying out and does not itself generate any particles. Wipes are available as 68 gsm binder-free, hydroentangled cellulose / polyester.

Product ordering

Format: Impregnated Pouch Wipes 420 x 250 mm

Case quantity: 10 x 20 wipes

Ecolab code: 6600554

[Product datasheet \(PDF\)](#)

Product Icon Attributes

Sterile, Low residue peroxide, Polyester / cellulose, 10 packs of 20 wipes

Brand

Ecolab

Type

Mop Wipe

Lead time (if out of stock)

15 days

Pack Size

200 (10x20)

Construction/Material

Non Woven

Size

42cmx25cm

Colour

White

Country of origin

United Kingdom

Sterility

Sterile



www.cleanroomshop.com
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com
Fax: + 44 (0) 1524 811589