

# Trespa Toplab Table (No Upstand and No Undershelf)



## Why Choose Trespa Toplab Table (No Upstand and No Undershelf)

The Trespa Toplab Base (Athlon) Tables are suitable for working environments where a high standard of hygiene and cleanliness is essential.

No upstand - centre table

No undershelf - leg room underneath

Trespa surface finish (colour - white E0-00) with mild steel powder coated legs (colour - white RAL 9010)

Ideal for use within cleanrooms

Extremely smooth surface ensuring reduced particle retention

Panels are easy to keep clean and as they do not absorb or retain moisture they can be cleaned frequently

No self assembly required

Made to order, bespoke sizes available, sturdy

## More Info

The Trespa Toplab Base surface structure is extremely smooth ensuring reduced particle retention

Panels are easy to keep clean and as they do not absorb or retain moisture they can be cleaned frequently

Made to order, bespoke sizes

The Trespa surface is sturdy, impact resistant, moisture resistant, scratch resistant and chemical resistant, as well as being antistatic

No self assembly required

Discounted prices available for orders of 6 or more tables

Manufactured from 16mm Trespa Toplab Base (Athlon) with 2mm chamfered edge finish with 20mm radius to outside corners

Includes heavy duty adjustable nylon feet - optional braked swivel non marking castors available



[www.cleanroomshop.com](http://www.cleanroomshop.com)

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)

Tel: + 44 (0) 1524 813022

Fax: + 44 (0) 1524 811589

30%

Fully packaged and protected for transportation purposes

Suitable for Cleanrooms, Laboratories, Aerospace, Defense, Electronics, Food Processing, Pharmaceutical Processing, Biotechnological Processing, Computer Assembly

Impact resistant, moisture resistant, scratch resistant and chemical resistant to the Trespa surface

Antistatic.

### Available Options

- 1200 x 600 x 900mm
- 1500 x 600 x 900mm
- 1500 x 700 x 900mm
- 1800 x 600 x 900mm
- 1800 x 700 x 900mm
- 1800 x 900 x 900mm

Available sizes	Weight
TT1 - 1200 x 600 x 900mm	31.05 kg
TT2 - 1500 x 600 x 900mm	36.63 kg
TT3 - 1500 x 700 x 900mm	40.70 kg
TT4 - 1800 x 600 x 900mm	42.21 kg
TT5 - 1800 x 700 x 900mm	46.94 kg
TT6 - 1800 x 900 x 900mm	56.40 kg

### Brand

Trespa

### Lead time (if out of stock)

8-10 weeks

### Pack Size

1

### Materials

Trespa - Cellulose Fibre 70%, Phenolic Resin 30%

### Construction/Material

Trespa surface

### Height

900mm

### Colour

Trespa - White (E0-00), Powder Coated Frame - White (RAL 9010)

### Applications

These tables are suitable for working environments where a high standard of hygiene and cleanliness is essential. Cleanrooms, Laboratories, Aerospace, Defense Electronics, Food Processing, Pharmaceutical Processing, Biotechnological Processing, Computer Assembly

### Technical data

Top is manufactured from 16mm Trespa Toplab Base with 2mm chamfered edge finish and 20mm radius to corners. Screwed to welded table frame manufactured from 40mm sq square section, powder coated. Trespa surface finish with

### Packaging

Fully packaged and protected for transportation purposes



www.cleanroomshop.com  
Tel: + 44 (0) 1524 813022

info@cleanroomshop.com  
Fax: + 44 (0) 1524 811589

## Sterility

Non-Sterile

## PPE Requirements

Precautions should be taken to protect the head, hands and feet

## Disposal

Incineration

## Table features

No Upstand and No Undershelf

## Loading capacity

150 kg/m

## Antistatic properties

1



[www.cleanroomshop.com](http://www.cleanroomshop.com)  
Tel: + 44 (0) 1524 813022

[info@cleanroomshop.com](mailto:info@cleanroomshop.com)  
Fax: + 44 (0) 1524 811589